**TECCare<sup>®</sup>**  
**CONTROL**

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## References

1. Eastmead L. Adopting a new cleaner disinfectant to optimise cleanliness and improve the staff experience. National Infection Prevention Society Conference, London 2013. Poster Presentation
2. Keward J. Disinfectants in health care: finding an alternative to chlorine dioxide. Br J Nurs. 2013 Sep 12-25; 22(16):926, 928-32.
3. Lowe S. Adopting a chlorine-free sporicidal wipe into clinical practice as part of a comprehensive care bundle to help reduce HCAI incidence. National Infection Prevention Society Conference, London 2013. Poster Presentation.
4. Embleton G. Improving process efficiency and medical device cleanliness at a foundation trust by adopting a new combined cleaner/disinfectant into routine use in a busy medical device library. National Infection Prevention Society Conference, Glasgow 2014. Poster Presentation.
5. Lees K. Using a new low temperature, non-corrosive laundry system to improve the cleanliness of high performance mattresses and reduce energy costs associated with the laundry process. National Infection Prevention Society Conference, Glasgow 2014. Poster Presentation.
6. Woodall C. Phase 3 testing on an ultra-high level sporicidal disinfectant: the missing link between laboratory based efficacy and clinical performance? National Infection Prevention Society Conference, Glasgow 2014. Poster Presentation.
7. Rutala WA, Weber DJ and the Healthcare Infection Control Practices Advisory Committee (HICPAC). Guideline for disinfection and sterilization in healthcare facilities, 2008. [homepage on the internet]; 2008. [cited 2013 Dec 10] Available from: [http://www.cdc.gov/hicpac/pdf/guidelines/disinfection\\_nov\\_2008.pdf](http://www.cdc.gov/hicpac/pdf/guidelines/disinfection_nov_2008.pdf).
8. Buck KM. Cleaning and disinfecting: The effects of germicides [homepage on the internet]. Infection Control Today. September 2001 [cited 2013 Dec 9]. Available from: <http://www.infectioncontrolday.com/articles/2001/09/infection-control-today-09-2001-the-effects-of-ge.aspx>.
9. Denyer SP, Stewart GSAB. Mechanisms of action of disinfectants. International Biodeterioration and Biodegradation 1998 Jan;41:261-8.
10. Gomi M, Osaki Y, Mori M, Sakagami Y. Synergistic bactericidal effects of a sublethal concentration of didecylmethylammonium chloride (DDAC) and low concentrations of nonionic surfactants against Staphylococcus aureus. Biocontrol Sci. 2012;17(4):175-81.
11. McDonnell G, Russell AD. Antiseptics and disinfectants: activity, action, and resistance. Clin Microbiol Rev. 1999 Jan;12(1):147-79.
12. Ioannou CJ, Hanlon GW, Denyer SP. Action of disinfectant quaternary ammonium compounds against Staphylococcus aureus. Antimicrob Agents Chemother. 2007 Jan;51(1):296-306.
13. Salton MR. Lytic agents, cell permeability, and monolayer penetrability. J Gen Physiol. 1968 Jul 1;52(1):227-52.
14. Maillard JY. Bacterial target sites for biocide action. J Appl Microbiol. 2002;92 Suppl:16S-27S.

15. Hospital Infection Research Laboratory. Efficacy tests (EN1276). EN1276 Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of chemical disinfectants and antiseptics used in food, industrial, domestic and institutional areas (phase 2/step 1). Hospital Infection Research Laboratory, Birmingham, UK. February 2008. Data on file, TECcare Antimicrobial Technologies.
16. Melbec Microbiology Laboratories (UKAS). EN efficacy test. Rossendale, Blackburn, UK. Data on file, TECcare Antimicrobial Technologies.
17. ATS Labs. AOAC Bactericidal efficacy test. Eagan, MN, USA. April 2008. Data on file, TECcare Antimicrobial Technologies.
18. BluTest Laboratories Ltd. Test Report: EN14476:2005 Chemical disinfectants and antiseptics – virucidal quantitative suspension test for chemical disinfectants and antiseptics used in human medicine - Test method and requirements (phase 2/step 1) Modified for Feline calicivirus (Human Norovirus Surrogate). Glasgow, UK. November 2012. Data on file, TECcare Antimicrobial Technologies.
19. Bioscience Laboratories Inc. AOAC Virucidal efficacy test. Bozeman, Montana (MA) USA. January 2009. Data on file, TECcare Antimicrobial Technologies.
20. Manchester University. BS EN 1650 Chemical disinfectants and antiseptics – quantitative suspension test for the evaluation of fungicidal activity of chemical disinfectants and antiseptics used in food, industrial, domestic and institutional areas. Manchester, UK. March 2007. Data on file, TECcare Antimicrobial Technologies.
21. ATS Labs. AOAC Fungicidal efficacy test. Eagan, MN, USA. April 2008. Data on file, TECcare Antimicrobial Technologies.
22. BluTest Laboratories Ltd. Test Report: EN13704:2002 Chemical disinfectants – quantitative suspension test for the evaluation of sporicidal activity of chemical disinfectants used in food, industrial, domestic and institutional areas – (phase 2, step 1), modified. Clostridium difficile endospores. Glasgow, UK. November 2012. Data on file, TECcare Antimicrobial Technologies.
23. BluTest Laboratories Ltd. Test Report: EN14348:2005 Chemical disinfectants and antiseptics – quantitative suspension test for the evaluation of mycobactericidal activity of chemical disinfectants in the medical area including instrument disinfectants – test methods and requirements (phase 2, step 1). Glasgow, UK. December 2014. Data on file, TECcare Antimicrobial Technologies.
24. Speight S, Moy A, Macken S, Chitnis R, Hoffman PN, Davies A, Bennett A, Walker JT. Evaluation of the sporicidal activity of different chemical disinfectants used in hospitals against Clostridium difficile. J Hosp Infect. 2011 Sep;79(1):18-22.

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**Kimmoisa Oy**
**TECcare Group Limited**